



January 2016

**Standard procedures**  
(abbreviated SOP's)

**used for the DNDi *in vitro* screening**

**against**

**Sleeping sickness**  
**Chagas disease**  
**Leishmaniasis**  
**Cytotoxicity**

**Prof. Dr. L. Maes, DVM, PhD**  
Laboratory of Microbiology, Parasitology and Hygiene (LMPH)  
University of Antwerp, campus CDE, S7.27  
Universiteitsplein 1  
B-2610 Wilrijk-Antwerpen  
Belgium

Tel: (32)3.265.26.27  
[louis.maes@ua.ac.be](mailto:louis.maes@ua.ac.be)



### Interpretation of results

- Activity is expressed as IC<sub>50</sub>-values in µM concentrations. For compounds without exact molecular weight, IC<sub>50</sub>-values are expressed in µg/ml.
- Based on the level of the IC<sub>50</sub>, semi-quantitative activity scores are given (normalized across the different assays to allow comparison)
- Compounds with score  $\geq 3$  may require confirmation testing and further follow-up **if** the observed activity is selective (absence of obvious cytotoxicity, high selectivity index).

### In vitro activity of reference compounds

UA	Compound code	Generic name	MW	Screen	IC <sub>50</sub> (µM)	
<b>R120</b>	UA/ K020	TAMOXIFEN CITRATE SALT	564	MRC-5 SV2	10.2	1.44
<b>R139</b>	TDR/ 42268/1	NICLOSAMIDE	327		7.11	3.71
<b>R126</b>		MILTEFOSIN	408	L.inf	10.07	4.23
<b>R116</b>	UA/ K016	AMPHOTERICIN B	924		1.16	1.25
<b>R126</b>		MILTEFOSIN	408	L.don	6.88	3.64
<b>R116</b>	UA/ K016	AMPHOTERICIN B	924		0.99	0.28
<b>R125</b>	TDR/ 10164/1	BENZNIDAZOL	260	T.cruz	2.89	0.56
<b>R132</b>	TDR/ 10739/1	NIFURTIMOX	287		0.91	0.41
<b>R131</b>	TDR/ 10738/1	SURAMIN	1297	T.b.bruc	0.03	0.01
<b>R152</b>	TDR/ 9957/1	MELARSOPROL	398		0.02	0.004
<b>R131</b>	TDR/ 10738/1	SURAMIN	1297	T.b.rhod	0.04	0.01
<b>R152</b>	TDR/ 9957/1	MELARSOPROL	398		0.005	0.002

*IC<sub>50</sub>-values of reference compounds are based on historical control values.*

*Only one reference drug is included per assay.*



<b>Toxicity</b>	<b><i>In vitro</i> cytotoxicity evaluation on human fibroblasts (MRC-5 cell line)</b>
-----------------	---

### Parasite and cell cultures

MRC-5<sub>SV2</sub><sup>1</sup> cells are cultured in MEM + Earl's salts-medium, supplemented with L-glutamine, NaHCO<sub>3</sub> and 5% inactivated fetal calf serum. All cultures and assays are conducted at 37°C under an atmosphere of 5% CO<sub>2</sub>. Other cell types (J774, L6, Vero, Hela, *e.a.*) can also be used for determination of cytotoxicity/selectivity.

### Compound solutions/dilutions

Compound stock solutions are prepared in 100% DMSO at 20 mM or mg/ml<sup>2</sup>. The compounds are serially pre-diluted (2-fold or 4-fold) in DMSO followed by a further (intermediate) dilution in demineralized water to assure a final in-test DMSO concentration of <1%.

### Drug sensitivity assays

Assays are performed in sterile 96-well microtiter plates, each well containing 10 µl of the watery compound dilutions together with 190 µl of MRC-5<sub>SV2</sub> inoculum (1.5x10<sup>5</sup> cells/ml). Cell growth is compared to untreated-control wells (100% cell growth) and medium-control wells (0% cell growth). After 3 days incubation, cell viability is assessed fluorimetrically after addition of 50µl resazurin per well<sup>3</sup>. After 4 hours at 37°C, fluorescence is measured ( $\lambda_{ex}$  550 nm,  $\lambda_{em}$  590 nm). The results are expressed as % reduction in cell growth/viability compared to control wells and an IC<sub>50</sub> and an IC<sub>90</sub> (50% and 90% inhibitory concentrations) are determined.

### Primary screen

The MRC-5<sub>SV2</sub> cell-line is used. The compounds are tested at 5 or 10 concentrations depending upon the number of compounds (4-fold compound dilutions) (64 - 16 - 4 - 1 - 0.25 - 0.0625 - 0.015625 - 0.0039 - 0.000975 and 0.00024 µM or µg/ml). The compound is classified non-toxic when the IC<sub>50</sub> is higher than 30 µM. Between 30 and 10 µM, the compound is regarded as moderately toxic. When the IC<sub>50</sub> is lower than 10 µM, the compound is classified as highly toxic. Cytotoxic reference compounds include vinblastine or paclitaxel (IC<sub>50</sub> <0.01µM), but these are rarely included because of health hazards for laboratory personnel. Other compounds are therefore used: tamoxifen or niclosamide.

### Secondary screen

The IC<sub>50</sub> is determined using an extended dose range (2-fold compound dilutions) still with a highest concentration of 64µM. Other cell lines and primary cells can be included: L6, J774, Hela, Vero and PMM (primary mouse macrophages). Parallel cytotoxicity evaluation is required to assess selective action of test compounds in parasite screens.

<sup>1</sup> MRC-5<sub>SV2</sub> cells are diploid human embryonic lung fibroblasts. A SV-40 transformed cell line is now available to obtain unlimited sub-cultivation characteristics.

<sup>2</sup> Concentrations are standard expressed in molar concentrations, except for natural products, drug mixtures and if the molecular weight is not known.

<sup>3</sup> Resazurin stock solution in phosphate buffer: 50µg/ml. Alamar Blue™ can be used as alternative: 5µl of a 1/10 Alamar Blue™ solution is added to each well



<b>Leishmania</b>	<b><i>In vitro</i> drug screening model against <i>Leishmania donovani</i> and <i>Leishmania infantum</i></b>
-------------------	---

### Parasite and cell cultures

Two *Leishmania* species (*L.infantum* MHOM/MA(BE)/67 and *L.donovani* MHOM/ET/67/L82) are used. The strains are maintained in the Golden Hamster (*Mesocricetus auratus*). Amastigotes are collected from the spleen of an infected donor hamster using three centrifugation purification steps (300 rpm, keeping the supernatants, 2200 rpm, keeping the supernatants and 3500 rpm, keeping the pellet) and spleen parasite burdens are assessed using the Stauber technique<sup>4</sup>. Primary peritoneal mouse macrophages are used as host cell and are collected 2 days after peritoneal stimulation with a 2% potato starch suspension. All cultures and assays are conducted at 37°C under an atmosphere of 5% CO<sub>2</sub>.

### Compound solutions/dilutions

Compound stock solutions are prepared in 100% DMSO at 20 mM or mg/ml<sup>5</sup>. The compounds are serially pre-diluted (2-fold or 4-fold) in DMSO followed by a further (intermediate) dilution in demineralized water to assure a final in-test DMSO concentration of <1%.

### Drug sensitivity assays

Assays are performed in 96-well microtiter plates, each well containing 10 µl of the compound dilutions together with 190 µl of macrophage/parasite inoculum (3.10<sup>4</sup> cells + 4.5.10<sup>5</sup> parasites/well). The inoculum is prepared in RPMI-1640 medium, supplemented with 200 mM L-glutamine, 16.5 mM NaHCO<sub>3</sub>, and 5% inactivated fetal calf serum. The macrophages are infected after 48 hours. The compounds are added after 2 hours of infection. Parasite multiplication is compared to untreated-infected controls (100% growth) and uninfected controls (0% growth). After 5 days incubation, parasite burdens (mean number of amastigotes/macrophage) are microscopically assessed after staining the cells with a 10% Giemsa solution. The results are expressed as % reduction in parasite burden compared to untreated control wells and an IC<sub>50</sub> and an IC<sub>90</sub> (50% and 90% inhibitory concentrations) are calculated.

### Primary screen

*Leishmania infantum* MHOM/MA(BE)/67 strain is used. The compounds are tested at 5 or 10 concentrations depending upon the number of compounds (4-fold compound dilutions) (64 – 16 – 4 – 1 – 0.25 – 0.0625 – 0.015625 – 0.0039 – 0.000975 and 0.00024 µM or µg/ml). Amphotericin B and miltefosin are included as the reference drugs. A test compound is classified as inactive when the IC<sub>50</sub> is higher than 30 µM. When IC<sub>50</sub> lies between 30 and 10 µM, the compound is regarded as moderately active. If the IC<sub>50</sub> is lower than 10 µM, the compound is classified as highly active on the condition that it also demonstrates selective action (absence of cytotoxicity against primary peritoneal macrophages). A final recommendation for activity is given after confirmatory evaluation in a secondary screening.

### Secondary screen

*Leishmania infantum* MHOM/MA(BE)/67 and *L.donovani* MHOM/ET/67/L82 strains are used and the IC<sub>50</sub>-values are determined using an extended dose range (2-fold compound dilutions). Amphotericin B and miltefosin are included as reference drugs.

<sup>4</sup> Stauber LA. (1966): Characterization of strains of *Leishmania donovani*. Exp Parasitol. 18: 1-11.

<sup>5</sup> Concentrations are standard expressed in molar concentrations, except for natural products, drug mixtures and if the molecular weight is not known.



<b>Chagas</b>	<b><i>In vitro</i> drug screening model against <i>Trypanosoma cruzi</i></b>
---------------	--

### **Parasite and cell cultures**

*Trypanosoma cruzi*, Tulahuen CL2,  $\beta$  galactosidase strain (nifurtimox-sensitive) is used<sup>6</sup>. The strain is maintained on MRC-5<sub>SV2</sub> (human lung fibroblast) cells<sup>7</sup> in MEM medium, supplemented with 200 mM L-glutamine, 16.5 mM NaHCO<sub>3</sub>, and 5% inactivated fetal calf serum. All cultures and assays are conducted at 37°C under an atmosphere of 5% CO<sub>2</sub>.

### **Compound solutions/dilutions**

Compound stock solutions are prepared in 100% DMSO at 20 mM or mg/ml<sup>8</sup>. The compounds are serially pre-diluted (2-fold or 4-fold) in DMSO followed by a further (intermediate) dilution in demineralized water to assure a final in-test DMSO concentration of <1%.

### **Drug sensitivity assays**

Assays are performed in sterile 96-well microtiter plates, each well containing 10  $\mu$ l of the watery compound dilutions together with 190  $\mu$ l of MRC-5 cell/parasite inoculum ( $4 \cdot 10^3$  cells/well +  $4 \cdot 10^4$  parasites/well). Parasite growth is compared to untreated-infected controls (100% growth) and non-infected controls (0% growth) after 7 days incubation at 37°C and 5% CO<sub>2</sub>. Parasite burdens are assessed after adding the substrate CPRG (chlorophenolred  $\beta$ -D-galactopyranoside): 50 $\mu$ l/well of a stock solution containing 15.2 mg CPRG + 250  $\mu$ l Nonidet in 100 ml PBS. The change in color is measured spectrophotometrically at 540 nm after 4 hours incubation at 37 °C. The results are expressed as % reduction in parasite burdens compared to control wells and an IC<sub>50</sub> and IC<sub>90</sub> (50% and 90% inhibitory concentrations) are calculated.

### **Primary screen**

*Trypanosoma cruzi*  $\beta$  galactosidase strain is used. Compounds are tested at 5 or 10 concentrations depending upon the number of compounds (4-fold compound dilutions) (64 – 16 – 4 – 1 – 0.25 – 0.0625 – 0.015625 – 0.0039 – 0.000975 and 0.00024  $\mu$ M or  $\mu$ g/ml). Nifurtimox or benznidazole can be included as the reference drugs. The test compound is classified as inactive when the IC<sub>50</sub> is higher than 30  $\mu$ M. When IC<sub>50</sub> lies between 30 and 5  $\mu$ M, the compound is regarded as being moderate active. When the IC<sub>50</sub> is lower than 5  $\mu$ M, the compound is classified as highly active on the condition that it also demonstrates selective action (absence of cytotoxicity). A final recommendation for activity is given after confirmatory evaluation in a secondary screening.

### **Secondary screen**

*Trypanosoma cruzi*  $\beta$  galactosidase strain is used and IC<sub>50</sub>-values are determined using an extended dose range (2-fold compound dilutions). Nifurtimox or benznidazole is included as reference drugs.

<sup>6</sup> Buckner FS, Verlinde CL, La Flamme AC, Van Voorhis WC. (1996): Efficient technique for screening drugs for activity against *Trypanosoma cruzi* using parasites expressing beta-galactosidase. Antimicrob Agents Chemother. 40: 2592-2597.

<sup>7</sup> MRC-5<sub>SV2</sub> cells are diploid human embryonic lung fibroblasts. A SV-40 transformed cell line is now available to obtain unlimited sub-cultivation characteristics.

<sup>8</sup> Concentrations are standard expressed in molar concentrations, except for natural products, drug mixtures and if the molecular weight is not known.



**Sleeping sickness**

***In vitro* drug screening model against *Trypanosoma brucei***

**Parasite and cell cultures**

The *Trypanosoma brucei brucei* Squib 427 strain (suramin-sensitive) or *Trypanosoma brucei rhodesiense* (strain STIB-900) are used for screening. The strains are maintained in Hirumi (HMI-9) medium, supplemented with 10% inactivated fetal calf serum. All cultures and assays are conducted at 37°C under an atmosphere of 5% CO<sub>2</sub>.

**Compound solutions/dilutions**

Compound stock solutions are prepared in 100% DMSO at 20 mM or mg/ml<sup>9</sup>. The compounds are serially pre-diluted (2-fold or 4-fold) in DMSO followed by a further (intermediate) dilution in demineralized water to assure a final in-test DMSO concentration of <1%.

**Drug sensitivity assays**

Assays are performed in sterile 96-well microtiter plates, each well containing 10 µl of the compound dilutions together with 190 µl of the parasite suspension (1.5x10<sup>4</sup> parasites/well - *T.brucei* or 4x10<sup>3</sup> parasites/well - *T.rhodesiense*). Parasite growth is compared to untreated-infected (100% parasite growth) and uninfected controls (0% growth). After 3 days incubation, parasite growth is assessed fluorimetrically after addition of 50µl resazurin per well<sup>10</sup>. After 6 hours (*T.rhodesiense*) or 24 hours (*T.brucei*) at 37°C, fluorescence is measured (λ<sub>ex</sub> 550 nm, λ<sub>em</sub> 590 nm). The results are expressed as % reduction in parasite growth/viability compared to control wells and an IC<sub>50</sub> and IC<sub>90</sub> (50% and 90% inhibitory concentrations) are calculated.

**Primary screen**

*Trypanosoma brucei brucei* Squib 427 strain and *Trypanosoma brucei rhodesiense* STIB-900 strain are used. Compounds are tested at 5 or 10 concentrations depending upon the number of compounds (4-fold compound dilutions) (64 – 16 – 4 – 1 – 0.25 – 0.0625 – 0.015625 – 0.0039 – 0.000975 and 0.00024 µM or µg/ml). Suramin or Melarsoprol are included as the reference drugs. The compound is classified as inactive when the IC<sub>50</sub> is higher than 5 µM. When IC<sub>50</sub> lies between 5 and 1 µM, the compound is regarded as being moderate active. When the IC<sub>50</sub> is lower than 1 µM, the compound is classified as highly active on the condition that it also demonstrates selective action (absence of cytotoxicity). A final recommendation for activity is given after confirmatory evaluation in a secondary screening.

**Secondary screen**

*Trypanosoma brucei brucei* and *T.b.rhodesiense* are both included used and IC<sub>50</sub>-values are determined using an extended dose range (2-fold compound dilutions). Suramin and melarsoprol are included as reference drugs.

<sup>9</sup> Concentrations are standard expressed in molar concentrations, except for natural products, drug mixtures and if the molecular weight is not known.

<sup>10</sup> Resazurin stock solution in phosphate buffer: 50µg/ml. // Alamar Blue™ can be used as alternative: 5µl of a 1/10 Alamar Blue™ solution is added to each well